

Onze infuuspompen



	Cadd Legacy One	Cadd Legacy PCA	Cadd Legacy Plus	Cadd Prizm Vip	Cadd Solis Vip	Asena GS	Crono 5 V2	Crono PCA 50	Crono S-Pid 50 V2	Graseby MS-26	T34L PCA	T34L standaard	Asena GW	Body-guard 323	Body-guard 323 Color Vision	Cadd MS-3	Paradigm	
Intermediair	cassette					spuit							zak			cartridge		
Toedieningswijze																		
Subcutaan		✓			✓		✓	✓	✓	✓	✓	✓					✓	✓
Perifeer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
CVC	✓		✓	✓	✓								✓	✓	✓			
PAC	✓		✓	✓	✓								✓	✓	✓			
PICC	✓		✓	✓	✓								✓	✓	✓			
Epiduraal						✓					✓	✓						
Spinaal					✓						✓	✓						
Type																		
Mobiele pomp	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓
Niet-mobiele pomp						✓							✓					
Doseringshoeveelheid																		
3 ml cartridge																	✓	✓
10 ml spuit						✓	✓			✓	✓	✓						
20 ml spuit						✓	✓			✓	✓	✓						
30 ml spuit						✓				✓	✓	✓						
50 ml spuit						✓		✓	✓		✓	✓						
50 ml cassette	✓	✓	✓	✓	✓													
100 ml cassette	✓	✓	✓	✓	✓													
250 ml cassette	✓	✓	✓	✓	✓													
Infuuszak													✓	✓	✓			
Fles													✓	✓	✓			
Container													✓	✓	✓			
Doseringsfrequentie																		
Bolus		✓			✓						✓		✓	✓	✓			
Continu ml/uur		✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓
Continu ml/24 uur	✓									✓								
Intermitterend			✓	✓	✓								✓	✓	✓			
PCA		✓		✓	✓			✓			✓			✓	✓			
Stappenprogramma					✓									✓	✓			
Veel gebruikt voor																		
Antibioticum	✓		✓	✓	✓								✓	✓	✓			
Bloedtransfusie*													✓	✓	✓			
Corticosteroïd			✓		✓								✓	✓	✓			
Cytostaticum*	✓		✓	✓	✓								✓	✓	✓			
Diureticum			✓		✓	✓								✓	✓			
Dobutamine			✓		✓									✓	✓			
Enzymtherapie	✓				✓								✓	✓	✓			
Flolan	✓				✓								✓	✓	✓			
Immunoglobuline								✓	✓				✓	✓	✓			
Pijnmedicatie		✓			✓		✓	✓		✓	✓	✓		✓	✓			
Remodulin IV	✓																	
Remodulin SC																✓	✓	
TPV*				✓	✓									✓	✓			
Veletri	✓				✓									✓	✓			
Vocht	✓		✓	✓	✓								✓	✓	✓			
Overige kenmerken																		
Bedieningsgemak	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓	✓	✓	✓	✓
Pompgewicht (g)	392	383	392	568	595	2.500	115	140	140	185	400	400	1.490	365	398	90	105	
Inclusief tas	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓
Verwisselbare batterij(en)	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓
Verwisselbare accu				✓	✓									✓	✓			
Niet-verwisselbare accu						✓					✓	✓	✓	✓	✓			
Netstroom	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓	✓			
Alarmpuncties																		
Batterij/accu leeg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Infusie gereed/pomp gestopt	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Lucht in de lijn	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Occlusiedruk	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reservoir leeg	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Systeemfout	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

* Hiervoor zijn speciale startpakketten beschikbaar.

Ons assortiment wordt regelmatig vernieuwd.